

CHAMPIONS for Podiatry



BIOLOGICS



TRILLIANT
SURGICAL



Bone Block

TALAR DOME

The **Talar Dome Bone Block** is derived from the superior talus, one of the strongest bones in the body. This unicortical graft's structure consists of a tightly-woven cancellous scaffold which is intended to provide increased strength and maximize bone remodeling potential.

- ▶ Processed for consistent and predictable sized blocks, available in medium and large
- ▶ Can be used for procedures such as Cotton or Evans osteotomies, as well as other foot or ankle reconstructive applications
- ▶ Freeze-dried, provided sterile via low dose, low temperature, gamma irradiation with a sterility assurance level (SAL) of 10^{-6}
- ▶ Packaged with a shelf life of 5 years and may be stored in ambient temperatures

Cancellous Matrix Comparison



Part Number	Size Dimension (L x W x D)
421041-TS	Medium, 14-18 x 14-18 x 10-14mm
421042-TS	Large, 20-24 x 20-24 x 15-20mm

For more information and detailed surgical techniques, visit www.trilliantsurgical.com

All biologic products are minimally manipulated and intended for homologous use only. Trilliant Surgical is a distributor for CellRight Technologies

www.trilliantsurgical.com | Phone: 800-495-2919 | Fax: 877-778-3864 | 727 North Shepherd Drive | Suite 100 | Houston, TX 77007