

### UDI Cross Reference Tool for Trilliant Surgical Products

This catalog contains applicable GTIN information for Trilliant Surgical products in order to assist in deriving required UDI information.

All product shipped from Trilliant Surgical in Houston, TX contains all the necessary UDI information on the respective product's package labeling; however, for non-sterile implants and instruments that are cleaned and sterilized prior to use, they are separated from their original label and packaging in order to undergo cleaning and sterilization wherein the original labeling may be discarded.

When product becomes separated from their original labeling, UDI information can be derived by combining available product identifying information, Product Identifiers (PI), on devices with a cross referencing tool containing associated Device Identifiers (DI).

*Note: Trilliant Surgical utilizes Lot Numbers as Product Identifiers.*

Within the GTIN catalog, for each product, there is an applicable GTIN, which is inclusive of the respective product's Device Identifier (DI), the portion of the UDI that identifies the specific version or model a device and the device labeler.

For example, the DI for part 101-00-001, (DI=00812926020075), can identify the part as a 3S Hemi Implant, Small, Non-Sterile and Trilliant Surgical as the device labeler for the part on global scale.

Together, the UDI information can be generated;  $DI + PI = UDI$

The format to generate a UDI (for Trilliant Surgical product) is as follows:  
**(01)<Device Identifier>(10)<Product Identifier>**

For example, if the part is 101-00-001 (DI=00812926020075) and the associated Lot is TSL001234, the UDI is:

**(01)00812926020075(10)TSL001234**

#### Important to Note:

- Copies of the most current reference catalog are available on Trilliant Surgical's website, <http://trilliantsurgical.com/regulatory/>
- Obtain the most recent version of this document before use to ensure information is current and accurate for each represented product
- This catalog may include active and inactive product; all products may not currently be offered and/or available at the time of use
- For any questions, please contact Trilliant Surgical Customer Service:

Phone: 800-495-2919

Email: [info@trilliantsurgical.com](mailto:info@trilliantsurgical.com)

Fax: 877-778-3864

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
101-00-001	3S Hemi Implant, Small, Non-Sterile	00812926020075
101-00-002	3S Hemi Implant, Medium, Non-Sterile	00812926020082
101-00-003	3S Hemi Implant, Large, Non-Sterile	00812926020099
101-00-004	3S Hemi Implant, XLarge, Non-Sterile	00812926020105
102-10-007	7mm Twist Subtalar Implant	00812926020235
102-10-008	8mm Twist Subtalar Implant	00812926020242
102-10-009	9mm Twist Subtalar Implant	00812926020259
102-10-010	10mm Twist Subtalar Implant	00812926020266
102-10-011	11mm Twist Subtalar Implant	00812926020273
102-10-012	12mm Twist Subtalar Implant	00812926020280
102-20-007	7mm Disco Subtalar Implant	00812926020297
102-20-008	8mm Disco Subtalar Implant	00812926020303
102-20-009	9mm Disco Subtalar Implant	00812926020310
102-20-010	10mm Disco Subtalar Implant	00812926020327
102-20-011	11mm Disco Subtalar Implant	00812926020334
102-20-012	12mm Disco Subtalar Implant	00812926020341
111-00-001	3S Hemi Non-Sterile Sterilization Tray	00812926020112
111-00-002	3S Hemi Handle w/ Impactor Cap	00812926020129
111-00-003	3S Hemi Sizer S/M, Non-Sterile	00812926020136
111-00-004	3S Hemi Sizer L/XL, Non-Sterile	00812926020143
111-00-005	3S Hemi, Non-Sterile, Awl	00812926020150
111-00-006	3S Hemi S/M Broach, Non-Sterile	00812926020167
111-00-007	3S Hemi L/XL Broach, Non-Sterile	00812926020174
111-00-008	3S Hemi Trial, Small, Non-Sterile	00812926020181
111-00-009	3S Hemi Trial, Medium, Non-Sterile	00812926020198
111-00-010	3S Hemi Trial, Large, Non-Sterile	00812926020204
111-00-011	3S Hemi Trial, XLarge, Non-Sterile	00812926020211
111-00-012	3S Hemi Impactor Bit, Non-Sterile	00812926020228
111-00-013	3S Hemi Impactor Head, Non-Sterile	00812926028033
111-99-002	3S Hemi, Non-Sterile, Tray Lid	00812926027852
112-00-001	Subtalar Implant System Sterilization Tray Base Assembly	00812926020358
112-00-002	Probe	00812926020365
112-00-003	Insertion Tool	00812926020372
112-00-004	Implant Pin	00812926020389
112-00-005	Removal Tool	00812926020396
112-00-006	Subtalar Tray Shipping Insert	00812926020402
112-10-007	7mm Twist Subtalar Trial	00812926020419
112-10-008	8mm Twist Subtalar Trial	00812926020426
112-10-009	9mm Twist Subtalar Trial	00812926020433
112-10-010	10mm Twist Subtalar Trial	00812926020440
112-10-011	11mm Twist Subtalar Trial	00812926020457
112-10-012	12mm Twist Subtalar Trial	00812926020464
112-20-007	7mm Disco Subtalar Trial	00812926020471
112-20-008	8mm Disco Subtalar Trial	00812926020488
112-20-009	9mm Disco Subtalar Trial	00812926020495
112-20-010	10mm Disco Subtalar Trial	00812926020501

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
112-20-011	11mm Disco Subtalar Trial	00812926020518
112-20-012	12mm Disco Subtalar Trial	00812926020525
112-99-011	Subtalar Implant System Sterilization Tray Lid	00812926027869
200-20-008	2.0mm x 8mm Tiger Cannulated Screw	00812926020532
200-20-010	2.0mm x 10mm Tiger Cannulated Screw	00812926020549
200-20-012	2.0mm x 12mm Tiger Cannulated Screw	00812926020556
200-20-014	2.0mm x 14mm Tiger Cannulated Screw	00812926020563
200-20-016	2.0mm x 16mm Tiger Cannulated Screw	00812926020570
200-20-018	2.0mm x 18mm Tiger Cannulated Screw	00812926020587
200-20-020	2.0mm x 20mm Tiger Cannulated Screw	00812926020594
200-20-022	2.0mm x 22mm Tiger Cannulated Screw	00812926020600
200-20-024	2.0mm x 24mm Tiger Cannulated Screw	00812926020617
200-20-026	2.0mm x 26mm Tiger Cannulated Screw	00812926020624
200-20-028	2.0mm x 28mm Tiger Cannulated Screw	00812926020631
200-20-030	2.0mm x 30mm Tiger Cannulated Screw	00812926020648
200-20-032	2.0mm x 32mm Tiger Cannulated Screw	00812926020655
200-20-034	2.0mm x 34mm Tiger Cannulated Screw	00812926020662
200-20-036	2.0mm x 36mm Tiger Cannulated Screw	00812926020679
200-20-038	2.0mm x 38mm Tiger Cannulated Screw	00812926020686
200-20-040	2.0mm x 40mm Tiger Cannulated Screw	00812926020693
200-20-042	2.0mm x 42mm Tiger Cannulated Screw	00812926020709
200-24-008	2.4mm x 8mm Tiger Cannulated Screw	00812926020716
200-24-010	2.4mm x 10mm Tiger Cannulated Screw	00812926020723
200-24-012	2.4mm x 12mm Tiger Cannulated Screw	00812926020730
200-24-014	2.4mm x 14mm Tiger Cannulated Screw	00812926020747
200-24-016	2.4mm x 16mm Tiger Cannulated Screw	00812926020754
200-24-018	2.4mm x 18mm Tiger Cannulated Screw	00812926020761
200-24-020	2.4mm x 20mm Tiger Cannulated Screw	00812926020778
200-24-022	2.4mm x 22mm Tiger Cannulated Screw	00812926020785
200-24-024	2.4mm x 24mm Tiger Cannulated Screw	00812926020792
200-24-026	2.4mm x 26mm Tiger Cannulated Screw	00812926020808
200-24-028	2.4mm x 28mm Tiger Cannulated Screw	00812926020815
200-24-030	2.4mm x 30mm Tiger Cannulated Screw	00812926020822
200-24-032	2.4mm x 32mm Tiger Cannulated Screw	00812926020839
200-24-034	2.4mm x 34mm Tiger Cannulated Screw	00812926020846
200-24-036	2.4mm x 36mm Tiger Cannulated Screw	00812926020853
200-24-038	2.4mm x 38mm Tiger Cannulated Screw	00812926020860
200-24-040	2.4mm x 40mm Tiger Cannulated Screw	00812926020877
200-24-042	2.4mm x 42mm Tiger Cannulated Screw	00812926020884
200-24-044	2.4mm x 44mm Tiger Cannulated Screw	00812926020891
200-24-046	2.4mm x 46mm Tiger Cannulated Screw	00812926020907
200-24-048	2.4mm x 48mm Tiger Cannulated Screw	00812926020914
200-30-010	3.0mm x 10mm Tiger Cannulated Screw	00812926020921
200-30-012	3.0mm x 12mm Tiger Cannulated Screw	00812926020938
200-30-014	3.0mm x 14mm Tiger Cannulated Screw	00812926020945
200-30-016	3.0mm x 16mm Tiger Cannulated Screw	00812926020952

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
200-30-018	3.0mm x 18mm Tiger Cannulated Screw	00812926020969
200-30-020	3.0mm x 20mm Tiger Cannulated Screw	00812926020976
200-30-022	3.0mm x 22mm Tiger Cannulated Screw	00812926020983
200-30-024	3.0mm x 24mm Tiger Cannulated Screw	00812926020990
200-30-026	3.0mm x 26mm Tiger Cannulated Screw	00812926021003
200-30-028	3.0mm x 28mm Tiger Cannulated Screw	00812926021010
200-30-030	3.0mm x 30mm Tiger Cannulated Screw	00812926021027
200-30-032	3.0mm x 32mm Tiger Cannulated Screw	00812926021034
200-30-034	3.0mm x 34mm Tiger Cannulated Screw	00812926021041
200-30-036	3.0mm x 36mm Tiger Cannulated Screw	00812926021058
200-30-038	3.0mm x 38mm Tiger Cannulated Screw	00812926021065
200-30-040	3.0mm x 40mm Tiger Cannulated Screw	00812926021072
200-40-010	4.0mm x 10mm Tiger Cannulated Screw	00812926021089
200-40-012	4.0mm x 12mm Tiger Cannulated Screw	00812926021096
200-40-014	4.0mm x 14mm Tiger Cannulated Screw	00812926021102
200-40-016	4.0mm x 16mm Tiger Cannulated Screw	00812926021119
200-40-018	4.0mm x 18mm Tiger Cannulated Screw	00812926021126
200-40-020	4.0mm x 20mm Tiger Cannulated Screw	00812926021133
200-40-022	4.0mm x 22mm Tiger Cannulated Screw	00812926021140
200-40-024	4.0mm x 24mm Tiger Cannulated Screw	00812926021157
200-40-026	4.0mm x 26mm Tiger Cannulated Screw	00812926021164
200-40-028	4.0mm x 28mm Tiger Cannulated Screw	00812926021171
200-40-030	4.0mm x 30mm Tiger Cannulated Screw	00812926021188
200-40-032	4.0mm x 32mm Tiger Cannulated Screw	00812926021195
200-40-034	4.0mm x 34mm Tiger Cannulated Screw	00812926021201
200-40-036	4.0mm x 36mm Tiger Cannulated Screw	00812926021218
200-40-038	4.0mm x 38mm Tiger Cannulated Screw	00812926021225
200-40-040	4.0mm x 40mm Tiger Cannulated Screw	00812926021232
200-40-042	4.0mm x 42mm Tiger Cannulated Screw	00812926021249
200-40-044	4.0mm x 44mm Tiger Cannulated Screw	00812926021256
200-40-046	4.0mm x 46mm Tiger Cannulated Screw	00812926021263
200-40-048	4.0mm x 48mm Tiger Cannulated Screw	00812926021270
200-40-050	4.0mm x 50mm Tiger Cannulated Screw	00812926021287
200-40-052	4.0mm x 52mm Tiger Cannulated Screw	00812926021294
200-40-054	4.0mm x 54mm Tiger Cannulated Screw	00812926021300
200-40-056	4.0mm x 56mm Tiger Cannulated Screw	00812926021317
200-45-001	4.5mm Smooth Washer	00812926029580
200-45-002	4.5mm Spiked Washer	00812926029573
200-45-020	4.5mm x 20mm STTCS	00812926028811
200-45-022	4.5mm x 22mm STTCS	00812926028828
200-45-024	4.5mm x 24mm STTCS	00812926028835
200-45-026	4.5mm x 26mm STTCS	00812926028842
200-45-028	4.5mm x 28mm STTCS	00812926028859
200-45-030	4.5mm x 30mm STTCS	00812926028866
200-45-032	4.5mm x 32mm STTCS	00812926028873
200-45-034	4.5mm x 34mm STTCS	00812926028880

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
200-45-036	4.5mm x 36mm STTCS	00812926028897
200-45-038	4.5mm x 38mm STTCS	00812926028903
200-45-040	4.5mm x 40mm STTCS	00812926028910
200-45-045	4.5mm x 45mm STTCS	00812926028927
200-45-050	4.5mm x 50mm STTCS	00812926028934
200-45-055	4.5mm x 55mm STTCS	00812926028941
200-45-060	4.5mm x 60mm STTCS	00812926028958
200-45-065	4.5mm x 65mm STTCS	00812926028965
200-45-070	4.5mm x 70mm STTCS	00812926028972
200-45-140	4.5mm x 40mm LTTCS	00812926028989
200-45-145	4.5mm x 45mm LTTCS	00812926028996
200-45-150	4.5mm x 50mm LTTCS	00812926029009
200-45-155	4.5mm x 55mm LTTCS	00812926029016
200-45-160	4.5mm x 60mm LTTCS	00812926029023
200-45-165	4.5mm x 65mm LTTCS	00812926029030
200-45-170	4.5mm x 70mm LTTCS	00812926029047
200-45-220	4.5mm x 20mm FTTCS	00812926029054
200-45-222	4.5mm x 22mm FTTCS	00812926029061
200-45-224	4.5mm x 24mm FTTCS	00812926029078
200-45-226	4.5mm x 26mm FTTCS	00812926029085
200-45-228	4.5mm x 28mm FTTCS	00812926029092
200-45-230	4.5mm x 30mm FTTCS	00812926029108
200-45-232	4.5mm x 32mm FTTCS	00812926029115
200-45-234	4.5mm x 34mm FTTCS	00812926029122
200-45-236	4.5mm x 36mm FTTCS	00812926029139
200-45-238	4.5mm x 38mm FTTCS	00812926029146
200-45-240	4.5mm x 40mm FTTCS	00812926029153
200-45-245	4.5mm x 45mm FTTCS	00812926029160
200-45-250	4.5mm x 50mm FTTCS	00812926029177
200-45-255	4.5mm x 55mm FTTCS	00812926029184
200-45-260	4.5mm x 60mm FTTCS	00812926029191
200-45-265	4.5mm x 65mm FTTCS	00812926029207
200-45-270	4.5mm x 70mm FTTCS	00812926029214
202-20-010	2.0 x 10mm Tiger Cann. Headless Screw	00812926021492
202-20-012	2.0 x 12mm Tiger Cann. Headless Screw	00812926021508
202-20-014	2.0 x 14mm Tiger Cann. Headless Screw	00812926021515
202-20-016	2.0 x 16mm Tiger Cann. Headless Screw	00812926021522
202-20-018	2.0 x 18mm Tiger Cann. Headless Screw	00812926021539
202-20-020	2.0 x 20mm Tiger Cann. Headless Screw	00812926021546
202-20-022	2.0 x 22mm Tiger Cann. Headless Screw	00812926021553
202-20-024	2.0 x 24mm Tiger Cann Headless Screw	00812926021560
202-24-010	2.4 x 10mm Tiger Cann. Headless Screw	00812926021577
202-24-012	2.4 x 12mm Tiger Cann. Headless Screw	00812926021584
202-24-014	2.4 x 14mm Tiger Cann. Headless Screw	00812926021591
202-24-016	2.4 x 16mm Tiger Cann. Headless Screw	00812926021607
202-24-018	2.4 x 18mm Tiger Cann. Headless Screw	00812926021614

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
202-24-020	2.4 x 20mm Tiger Cann. Headless Screw	00812926021621
202-24-022	2.4 x 22mm Tiger Cann. Headless Screw	00812926021638
202-24-024	2.4 x 24mm Tiger Cann. Headless Screw	00812926021645
202-24-026	2.4 x 26mm Tiger Cann. Headless Screw	00812926021652
202-24-028	2.4 x 28mm Tiger Cann. Headless Screw	00812926021669
202-24-030	2.4 x 30mm Tiger Cann. Headless Screw	00812926021676
202-30-010	3.0 x 10mm Tiger Cann. Headless Screw	00812926021683
202-30-012	3.0 x 12mm Tiger Cann. Headless Screw	00812926021690
202-30-014	3.0 x 14mm Tiger Cann. Headless Screw	00812926021706
202-30-016	3.0 x 16mm Tiger Cann. Headless Screw	00812926021713
202-30-018	3.0 x 18mm Tiger Cann. Headless Screw	00812926021720
202-30-020	3.0 x 20mm Tiger Cann. Headless Screw	00812926021737
202-30-022	3.0 x 22mm Tiger Cann. Headless Screw	00812926021744
202-30-024	3.0 x 24mm Tiger Cann. Headless Screw	00812926021751
202-30-026	3.0 x 26mm Tiger Cann. Headless Screw	00812926021768
202-30-028	3.0 x 28mm Tiger Cann. Headless Screw	00812926021775
202-30-030	3.0 x 30mm Tiger Cann. Headless Screw	00812926021782
202-30-032	3.0 x 32mm Tiger Cann. Headless Screw	00812926021799
202-30-034	3.0 x 34mm Tiger Cann. Headless Screw	00812926021805
202-45-020	4.5mm x 20mm Headless, STTHCS	00812926029221
202-45-022	4.5mm x 22mm Headless, STTHCS	00812926029238
202-45-024	4.5mm x 24mm Headless, STTHCS	00812926029245
202-45-026	4.5mm x 26mm Headless, STTHCS	00812926029252
202-45-028	4.5mm x 28mm Headless, STTHCS	00812926029269
202-45-030	4.5mm x 30mm Headless, STTHCS	00812926029276
202-45-032	4.5mm x 32mm Headless, STTHCS	00812926029283
202-45-034	4.5mm x 34mm Headless, STTHCS	00812926029290
202-45-036	4.5mm x 36mm Headless, STTHCS	00812926029306
202-45-038	4.5mm x 38mm Headless, STTHCS	00812926029313
202-45-040	4.5mm x 40mm Headless, STTHCS	00812926029320
202-45-045	4.5mm x 45mm Headless, STTHCS	00812926029337
202-45-050	4.5mm x 50mm Headless, STTHCS	00812926029344
202-45-055	4.5mm x 55mm Headless, STTHCS	00812926029351
202-45-060	4.5mm x 60mm Headless, STTHCS	00812926029368
202-45-065	4.5mm x 65mm Headless, STTHCS	00812926029375
202-45-070	4.5mm x 70mm Headless, STTHCS	00812926029382
204-24-006	2.6 x 2.4 x 6mm Hammertoe Implant	00812926028699
204-24-008	2.6 x 2.4 x 8mm Hammertoe Implant	00812926028705
204-24-010	2.6 x 2.4 x 10mm Hammertoe Implant	00812926028712
204-24-012	2.6 x 2.4 x 12mm Hammertoe Implant	00812926028729
204-24-014	2.6 x 2.4 x 14mm Hammertoe Implant	00812926028736
204-24-016	2.6 x 2.4 x 16mm Hammertoe Implant	00812926028743
204-24-018	2.6 x 2.4 x 18mm Hammertoe Implant	00812926028750
204-30-006	3.4 x 3.0 x 6mm Hammertoe Implant	00812926021935
204-30-008	3.4 x 3.0 x 8mm Hammertoe Implant	00812926021942
204-30-010	3.4 x 3.0 x 10mm Hammertoe Implant	00812926021959

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
204-30-012	3.4 x 3.0 x 12mm Hammertoe Implant	00812926021966
204-30-014	3.4 x 3.0 x 14mm Hammertoe Implant	00812926021973
204-30-016	3.4 x 3.0 x 16mm Hammertoe Implant	00812926021980
204-30-018	3.4 x 3.0 x 18mm Hammertoe Implant	00812926021997
204-30-020	3.4 x 3.0 x 20mm Hammertoe Implant	00812926022000
204-30-022	3.4 x 3.0 x 22mm Hammertoe Implant	00812926022017
204-30-024	3.4 x 3.0 x 24mm Hammertoe Implant	00812926022024
204-30-026	3.4 x 3.0 x 26mm Hammertoe Implant	00812926022031
204-30-028	3.4 x 3.0 x 28mm Hammertoe Implant	00812926022048
204-30-030	3.4 x 3.0 x 30mm Hammertoe Implant	00812926022055
204-30-106	4.5 x 3.0 x 6mm Hammertoe Implant	00812926022062
204-30-108	4.5 x 3.0 x 8mm Hammertoe Implant	00812926022079
204-30-110	4.5 x 3.0 x 10mm Hammertoe Implant	00812926022086
204-30-112	4.5 x 3.0 x 12mm Hammertoe Implant	00812926022093
204-30-114	4.5 x 3.0 x 14mm Hammertoe Implant	00812926022109
204-30-116	4.5 x 3.0 x 16mm Hammertoe Implant	00812926022116
204-30-118	4.5 x 3.0 x 18mm Hammertoe Implant	00812926022123
204-30-120	4.5 x 3.0 x 20mm Hammertoe Implant	00812926022130
204-30-122	4.5 x 3.0 x 22mm Hammertoe Implant	00812926022147
204-30-124	4.5 x 3.0 x 24mm Hammertoe Implant	00812926022154
204-30-126	4.5 x 3.0 x 26mm Hammertoe Implant	00812926022161
204-30-128	4.5 x 3.0 x 28mm Hammertoe Implant	00812926022178
204-30-130	4.5 x 3.0 x 30mm Hammertoe Implant	00812926022185
204-40-206	5.5mm x 4.0mm x 6mm Hammertoe Implant	00812926022192
204-40-208	5.5mm x 4.0mm x 8mm Hammertoe Implant	00812926022208
204-40-210	5.5mm x 4.0 x 10mm Hammertoe Implant	00812926022215
205-55-001	5.5mm Smooth Washer	00812926025070
205-55-030	5.5mm x 30mm STTCS	00812926024073
205-55-035	5.5mm x 35mm STTCS	00812926024080
205-55-040	5.5mm x 40mm STTCS	00812926024097
205-55-045	5.5mm x 45mm STTCS	00812926024103
205-55-050	5.5mm x 50mm STTCS	00812926024110
205-55-055	5.5mm x 55mm STTCS	00812926024127
205-55-060	5.5mm x 60mm STTCS	00812926024134
205-55-065	5.5mm x 65mm STTCS	00812926024141
205-55-070	5.5mm x 70mm STTCS	00812926024158
205-55-075	5.5mm x 75mm STTCS	00812926024165
205-55-080	5.5mm x 80mm STTCS	00812926024172
205-70-001	7.0mm Smooth Washer	00812926025087
205-70-050	7.0mm x 50mm STTCS	00812926024189
205-70-055	7.0mm x 55mm STTCS	00812926024196
205-70-060	7.0mm x 60mm STTCS	00812926024202
205-70-065	7.0mm x 65mm STTCS	00812926024219
205-70-070	7.0mm x 70mm STTCS	00812926024226
205-70-075	7.0mm x 75mm STTCS	00812926024233
205-70-080	7.0mm x 80mm STTCS	00812926024240

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
205-70-085	7.0mm x 85mm STTCS	00812926024257
205-70-090	7.0mm x 90mm STTCS	00812926024264
205-70-095	7.0mm x 95mm STTCS	00812926024271
205-70-100	7.0mm x 100mm STTCS	00812926024288
205-70-105	7.0mm x 105mm STTCS	00812926024295
205-70-110	7.0mm x 110mm STTCS	00812926024301
205-70-115	7.0mm x 115mm STTCS	00812926024318
205-70-120	7.0mm x 120mm STTCS	00812926024325
206-55-050	5.5mm x 50mm LTTCS	00812926024332
206-55-055	5.5mm x 55mm LTTCS	00812926024349
206-55-060	5.5mm x 60mm LTTCS	00812926024356
206-55-065	5.5mm x 65mm LTTCS	00812926024363
206-55-070	5.5mm x 70mm LTTCS	00812926024370
206-55-075	5.5mm x 75mm LTTCS	00812926024387
206-55-080	5.5mm x 80mm LTTCS	00812926024394
206-55-085	5.5mm x 85mm LTTCS	00812926024400
206-55-090	5.5mm x 90mm LTTCS	00812926024417
206-70-050	7.0mm x 50mm LTTCS	00812926024424
206-70-055	7.0mm x 55mm LTTCS	00812926024431
206-70-060	7.0mm x 60mm LTTCS	00812926024448
206-70-065	7.0mm x 65mm LTTCS	00812926024455
206-70-070	7.0mm x 70mm LTTCS	00812926024462
206-70-075	7.0mm x 75mm LTTCS	00812926024479
206-70-080	7.0mm x 80mm LTTCS	00812926024486
206-70-085	7.0mm x 85mm LTTCS	00812926024493
206-70-090	7.0mm x 90mm LTTCS	00812926024509
206-70-095	7.0mm x 95mm LTTCS	00812926024516
206-70-100	7.0mm x 100mm LTTCS	00812926024523
206-70-105	7.0mm x 105mm LTTCS	00812926024530
206-70-110	7.0mm x 110mm LTTCS	00812926024547
206-70-115	7.0mm x 115mm LTTCS	00812926024554
206-70-120	7.0mm x 120mm LTTCS	00812926024561
207-55-030	5.5 x 30mm STTHCS	00812926024578
207-55-035	5.5 x 35mm STTHCS	00812926024585
207-55-040	5.5 x 40mm STTHCS	00812926024592
207-55-045	5.5 x 45mm STTHCS	00812926024608
207-55-050	5.5 x 50mm STTHCS	00812926024615
207-55-055	5.5 x 55mm STTHCS	00812926024622
207-55-060	5.5 x 60mm STTHCS	00812926024639
207-55-065	5.5 x 65mm STTHCS	00812926024646
207-55-070	5.5 x 70mm STTHCS	00812926024653
207-55-075	5.5 x 75mm STTHCS	00812926024660
207-55-080	5.5 x 80mm STTHCS	00812926024677
207-70-045	7.0 x 45mm STTHCS	00812926027494
207-70-050	7.0 x 50mm STTHCS	00812926024684
207-70-055	7.0 x 55mm STTHCS	00812926024691



## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
207-70-060	7.0 x 60mm STTHCS	00812926024707
207-70-065	7.0 x 65mm STTHCS	00812926024714
207-70-070	7.0 x 70mm STTHCS	00812926024721
207-70-075	7.0 x 75mm STTHCS	00812926024738
207-70-080	7.0 x 80mm STTHCS	00812926024745
207-70-085	7.0 x 85mm STTHCS	00812926024752
207-70-090	7.0 x 90mm STTHCS	00812926024769
207-70-095	7.0 x 95mm STTHCS	00812926024776
207-70-100	7.0 x 100mm STTHCS	00812926024783
207-70-105	7.0 x 105mm STTHCS	00812926024790
207-70-110	7.0 x 110mm STTHCS	00812926024806
207-70-115	7.0 x 115mm STTHCS	00812926024813
207-70-120	7.0 x 120mm STTHCS	00812926024820
208-55-050	5.5 x 50mm LTHCS	00812926024837
208-55-055	5.5 x 55mm LTHCS	00812926024844
208-55-060	5.5 x 60mm LTHCS	00812926024851
208-55-065	5.5 x 65mm LTHCS	00812926024868
208-55-070	5.5 x 70mm LTHCS	00812926024875
208-55-075	5.5 x 75mm LTHCS	00812926024882
208-55-080	5.5 x 80mm LTHCS	00812926024899
208-55-085	5.5 x 85mm LTHCS	00812926024905
208-55-090	5.5 x 90mm LTHCS	00812926024912
208-70-050	7.0 x 50mm LTHCS	00812926024929
208-70-055	7.0 x 55mm LTHCS	00812926024936
208-70-060	7.0 x 60mm LTHCS	00812926024943
208-70-065	7.0 x 65mm LTHCS	00812926024950
208-70-070	7.0 x 70mm LTHCS	00812926024967
208-70-075	7.0 x 75mm LTHCS	00812926024974
208-70-080	7.0 x 80mm LTHCS	00812926024981
208-70-085	7.0 x 85mm LTHCS	00812926024998
208-70-090	7.0 x 90mm LTHCS	00812926025001
208-70-095	7.0 x 95mm LTHCS	00812926025018
208-70-100	7.0 x 100mm LTHCS	00812926025025
208-70-105	7.0 x 105mm LTHCS	00812926025032
208-70-110	7.0 x 110mm LTHCS	00812926025049
208-70-115	7.0 x 115mm LTHCS	00812926025056
208-70-120	7.0 x 120mm LTHCS	00812926025063
210-00-001	Tiger Screw Sterilization Tray Base Assembly	00812926021324
210-00-003	Cannulated Driver Handle	00812926021331
210-00-004	Ratcheting Cannulated Driver Handle	00812926021348
210-00-005	Cannulated Depth Gauge	00812926021355
210-00-006	Screw Remover	00812926021362
210-00-007	Osteotomy Guide	00812926027012
210-00-008	4.5mm Screw Remover	00812926029504
210-00-101	Combined Tiger / Tiger Headless Screw Tray	00812926027029
210-00-102	Compact Tiger Cannulated Screw Tray Base Assembly	00812926027036

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
210-00-103	4.5mm Tiger Cannulated Screw System Tray	00812926029535
210-24-001	2.0/2.4mm Cannulated Drill Bit	00812926021379
210-24-002	2.0/2.4mm Cannulated Countersink	00812926021386
210-24-003	2.0/2.4mm Cannulated Screw Driver Bit	00812926021393
210-24-004	0.035" x 5" K-Wire Standard	00812926021409
210-24-005	0.035" x 5" K-Wire Partially Threaded	00812926021416
210-24-006	0.035" x 5" K-Wire Double Ended	00812926021423
210-24-007	2.0/2.4mm Tiger Solid Core Slip Fit Driver	00812926029399
210-24-008	2.0/2.4mm Tiger Cannulated Slip Fit Driver	00812926029405
210-24-014	.035" Marked K-wire	00812926027821
210-30-002	Combo Depth Gauge Countersink	00812926027814
210-40-001	3.0/4.0mm Cannulated Drill Bit	00812926021430
210-40-002	3.0/4.0mm Cann. Screw Countersink	00812926021447
210-40-003	3.0/4.0mm Cannulated Screw Driver Bit	00812926021454
210-40-004	0.045" x 6" K-Wire Standard	00812926021461
210-40-005	0.045" x 6" K-Wire Partially Threaded	00812926021478
210-40-006	0.045" x 6" K-Wire Double Ended	00812926021485
210-40-007	3.0/4.0mm Tiger Solid Core Slip Fit Driver	00812926029412
210-40-008	3.0/4.0mm Tiger Cannulated Slip Fit Driver	00812926029429
210-40-014	0.045" x 6" Marked K-Wire Standard	00812926027838
210-45-001	4.5mm Cannulated Drill Bit	00812926029436
210-45-002	4.5mm Cannulated Countersink	00812926029450
210-45-003	4.5mm Cannulated Driver	00812926029467
210-45-005	4.5mm Cannulated Depth Gauge	00812926029498
210-45-006	4.5mm Drill Sleeve Assembly	00812926029481
210-45-007	4.5mm Solid-Core Screw Remover	00812926029511
210-45-011	4.5mm Cannulated Over Drill Bit	00812926029443
210-54-006	0.054" x 6" K-Wire Double Ended	00812926022314
210-60-004	0.062" x 6" K-Wire Standard	00812926023946
210-60-005	0.062" x 6" K-Wire Partially Threaded	00812926023953
210-60-006	0.062" x 6" K-Wire Double Ended	00812926023960
210-99-011	Tiger Cannulated Screw Tray Lid	00812926027876
210-99-101	Combined Tiger / Tiger Headless Screw Tray Lid	00812926027883
210-99-103	0.062" K-Wire Tube	00810028391369
210-99-112	Tiger Screw Compact Lid	00812926027890
212-00-001	Tiger Headless Screw Sterilization Tray Base Assembly	00812926021812
212-00-005	Cannulated Depth Gauge: Headless	00812926021829
212-20-002	2.0mm Cann. Prox. Drill	00812926021836
212-20-003	2.0mm Cann. Screw Driver Bit: Headless	00812926021843
212-20-007	2.0mm Removal Bit: Headless	00812926021850
212-24-002	2.4mm Cann. Prox. Drill	00812926021867
212-24-003	2.4mm Cann. Screw Driver Bit: Headless	00812926021874
212-24-007	2.4mm Removal Bit: Headless	00812926021881
212-30-002	3.0mm Cann. Prox. Drill	00812926021898
212-40-003	3.0/4.0mm Cann. Screw Driver Bit: Headless	00812926021911
212-40-007	3.0mm Removal Bit: Headless	00812926021928

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
212-45-002	4.5mm Cannulated Proximal Drill	00812926029474
212-99-011	Tiger Headless Cannulated Screw Tray Lid	00812926027906
214-00-001	Hammer Toe Implant Tray Base Assembly	00812926022222
214-00-006	3.4mm Spade Pilot Drill	00812926022239
214-00-007	4.5mm Spade Pilot Drill	00812926022246
214-00-008	5.5mm Spade Pilot Drill	00812926022253
214-00-009	3.4mm Spade Driver	00812926022260
214-00-010	5.5mm Spade Driver	00812926022277
214-00-011	4.5mm Spade Driver	00812926022284
214-00-012	Two-Step Cannulated Depth Gauge	00812926022291
214-00-013	Two-Step 3.0mm Thread Pilot Drill	00812926022307
214-00-014	Two-Step Hammer Toe Impactor	00812926027043
214-00-015	2.6mm Spade Pilot Drill	00812926028767
214-00-016	2.6mm Spade Driver	00812926028774
214-00-017	Two-Step 2.4mm Thread Pilot Drill	00812926028781
214-99-004	Hammer Toe Implant Tray Lid	00812926027913
214-99-008	2.6mm x 2.4mm Two Step Pocket Insert	00812926029931
215-00-001	Tiger Large Cannulated Screw Tray Assembly	00812926025094
215-00-003	Large Cannulated Fixed Driver Handle	00812926025117
215-00-004	Large Cannulated Ratcheting Driver Handle	00812926025124
215-00-005	Large Cannulated Depth Gauge	00812926025131
215-00-006	Large Cannulated Screw Sleeve	00812926025148
215-00-010	Large Cannulated Parallel Wire Guide	00812926025216
215-00-011	Large Cannulated Screw Extractor	00812926025223
215-00-012	1/4 SQ / Jacobs Adapter	00812926027050
215-55-001	5.5mm Large Cannulated Pilot Drill	00812926025230
215-55-002	5.5mm Large Cann. Countersink	00812926025247
215-55-003	5.5mm Large Cannulated Driver	00812926025254
215-55-004	0.094" x 12" K-Wire Standard	00812926027067
215-55-005	0.094" x 12" K-Wire Partially Threaded	00812926025315
215-55-007	5.5mm Screw Drill Sleeve	00812926025155
215-55-008	5.5mm Guide Pin Sleeve	00812926025162
215-55-009	5.5mm Screw Trocar	00812926025179
215-70-001	7.0mm Large Cannulated Pilot Drill	00812926025278
215-70-002	7.0mm Large Cann. Countersink	00812926025285
215-70-003	7.0mm Large Cannulated Driver	00812926025292
215-70-004	0.110" x 12" K-Wire Standard	00812926027074
215-70-005	0.110" x 12" K-Wire Partially Threaded	00812926025322
215-70-007	7.0mm Screw Drill Sleeve	00812926025186
215-70-008	7.0mm Guide Pin Sleeve	00812926025193
215-70-009	7.0mm Screw Trocar	00812926025209
217-55-002	5.5mm Cannulated Proximal Drill	00812926025261
217-55-003	5.5mm Headless Cannulated Driver	00812926027470
217-70-002	7.0mm Cannulated Proximal Drill	00812926025308
217-70-003	7.0mm Headless Cannulated Driver	00812926027487
220-00-001	Reamer Sterilization Tray	00812926027081

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
220-08-001	8mm Concave Reamer	00812926022321
220-08-002	8mm Convex Reamer	00812926022338
220-08-003	8mm Sterile HTR Kit	00812926022345
220-08-004	8mm HTR Implant System, Non-Sterile	00812926022352
220-10-001	10mm CONCAVE REAMER	00812926022376
220-10-002	10mm CONVEX REAMER	00812926022383
220-10-003	10mm Sterile HTR Kit	00812926022390
220-10-004	HTR Implant, Non-Sterile	00812926022406
220-12-001	12mm Concave Reamer	00812926027098
220-12-002	12mm Convex Reamer	00812926027104
220-12-004	12mm HTR Implant System, Non-Sterile	00812926027128
300-50-001	2 Hole VLC Gridlock Plate	00812926022420
300-50-002	3 Hole VLC Gridlock Plate	00812926022437
300-50-003	4 Hole VLC Gridlock Plate	00812926022444
300-50-004	Custom 4 Hole VLC Gridlock Plate	00812926022451
300-50-005	5 Hole VLC Gridlock Plate	00812926022468
300-50-006	6 Hole VLC Gridlock Plate	00812926022475
300-50-007	7 Hole VLC Gridlock Plate	00812926025339
300-50-008	Custom 5 Hole VLC Gridlock Plate	00812926025346
300-51-001	Small MPX Gridlock Plate	00812926022482
300-51-002	Medium MPX Gridlock Plate	00812926022499
300-51-003	Large MPX Gridlock Plate	00812926022505
300-52-001	Stepped MPJ VLC Gridlock Plate, Right	00812926022512
300-52-002	Stepped MPJ VLC Gridlock Plate, Left	00812926022529
300-52-003	Medium MPJ VLC Gridlock Plate	00812926022536
300-52-004	Petite MPJ VLC Gridlock Plate Right	00812926025353
300-52-005	Petite MPJ VLC Gridlock Plate Left	00812926025360
300-53-001	Lapidus VLC Gridlock Plate Right	00812926022543
300-53-002	Lapidus VLC Gridlock Plate Left	00812926022550
300-53-003	Lapidus VLC Gridlock Plate	00812926022567
300-54-002	2mm Wedge Gridlock Plate	00812926022574
300-54-004	4mm Wedge Gridlock Plate	00812926022581
300-54-006	6mm Wedge Gridlock Plate	00812926022598
300-54-102	2mm Pre-Bent Wedge VL Gridlock Plate	00812926022604
300-54-104	4mm Pre-Bent Wedge VL Gridlock Plate	00812926022611
300-54-106	6mm Pre-Bent Wedge VL Gridlock Plate	00812926022628
300-54-208	8mm Evans Plate	00812926027197
300-54-210	10mm Evans Plate	00812926027203
300-54-212	12mm Evans Plate	00812926027210
300-55-003	4 Hole VL Gridlock Plate	00812926027227
300-56-002	Small VLC "L" Plate Right	00812926025377
300-56-003	Small VLC "L" Plate Left	00812926025384
300-56-004	Large VLC "L" Plate Right	00812926025391
300-56-005	Large VLC "L" Plate Left	00812926025407
300-60-001	7 Hole VL Gridlock Fibula Plate	00812926025438
300-60-002	9 Hole VL Gridlock Fibula Plate	00812926025445

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
300-60-003	11 Hole VL Gridlock Fibula Plate	00812926025452
300-61-001	7 Hole Contoured VL Gridlock Fibula Plate	00812926025469
300-61-002	9 Hole Contoured VL Gridlock Fibula Plate	00812926025476
300-61-003	11 Hole Contoured VL Gridlock Fibula Plate	00812926025483
300-61-004	7 Hole VL Gridlock Distal Fibula Plate	00812926027562
300-61-005	8 Hole Contoured VL Gridlock Fibula Plate	00812926027579
300-61-006	10 Hole Contoured VL Gridlock Fibula Plate	00812926027586
300-61-007	12 Hole Contoured Fibula Plate	00812926027593
300-62-001	7 Hole Contour Tibia Plate	00812926025490
300-62-002	9 Hole Contour Tibia Plate	00812926025506
300-63-001	6 Hole 1/3 Tubular Plate	00812926027609
300-63-002	9 Hole 1/3 Tubular Plate	00812926027616
300-63-003	12 Hole 1/3 Tubular Plate	00812926027623
300-64-001	6 Hole Malleolar Gridlock Hook Plate	00812926027647
300-64-002	Buttress Plate	00812926027630
300-65-001	8 Hole Gridlock Helical Plate, Left	00812926028484
300-65-002	10 Hole Gridlock Helical Plate, Left	00812926028491
300-65-003	8 Hole Gridlock Helical Plate, Right	00812926028507
300-65-004	10 Hole Gridlock Helical Plate, Right	00812926028514
300-70-001	Minimally Invasive Bunion Plate, Right	00812926027708
300-70-002	Minimally Invasive Bunion Plate, Left	00812926027715
300-80-001	Arsenal Custom 4 Hole VLC Plate	00812926029948
300-80-002	Arsenal 2 Hole VLC Plate	00812926029955
300-80-003	Arsenal 3 Hole VLC Plate	00812926029962
300-80-004	Arsenal 4 Hole VLC Plate	00812926029979
300-80-005	Arsenal 5 Hole VLC Plate	00812926029986
300-80-006	Arsenal 6 Hole VLC Plate	00812926029993
300-81-001	Arsenal Mini Condylar Plate, Short	00810028390003
300-81-002	Arsenal Mini Condylar Plate, Long	00810028390010
300-81-003	Arsenal 3 Hole Mini VLC Plate	00810028390027
300-81-004	Arsenal 4 Hole Mini VLC Plate	00810028390034
300-81-005	Arsenal 5 Hole Mini VLC Plate	00810028390041
300-81-006	Arsenal 6 Hole Mini VLC Plate	00810028390058
300-81-007	Arsenal Mini Z Plate, Short	00810028390065
300-81-008	Arsenal Mini Z Plate, Long	00810028390072
300-82-001	Arsenal Petite MPJ Plate, L	00810028390089
300-82-002	Arsenal Petite MPJ Plate, R	00810028390096
300-82-003	Arsenal Standard MPJ Plate, L	00810028390102
300-82-004	Arsenal Standard MPJ Plate, R	00810028390119
300-82-005	Arsenal Revision MPJ Plate, L	00810028390126
300-82-006	Arsenal Revision MPJ Plate, R	00810028390133
300-83-001	Arsenal Petite Lapidus Plate, L	00810028390140
300-83-002	Arsenal Petite Lapidus Plate, R	00810028390157
300-83-003	Arsenal Standard Lapidus Plate, L	00810028390164
300-83-004	Arsenal Standard Lapidus Plate, R	00810028390171
300-83-005	Arsenal Large Lapidus Plate, L	00810028390188

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
300-83-006	Arsenal Large Lapidus Plate, R	00810028390195
300-84-001	Arsenal 5th Met Hook Plate, L	00810028390201
300-84-002	Arsenal 5th Met Hook Plate, R	00810028390218
300-85-001	Arsenal Lisfranc Plate	00810028390225
300-85-002	Arsenal TN Plate	00810028390232
300-85-003	Arsenal NC Plate, Left	00810028390249
300-85-004	Arsenal NC Plate, Right	00810028391345
300-86-001	Arsenal Evans Plate	00810028390256
300-86-002	Arsenal Universal Calc Slide Plate	00810028390263
300-86-003	Arsenal Calc Perimeter Plate, Left	00812926029825
300-86-004	Arsenal Calc Perimeter Plate, Right	00812926029832
300-86-005	Arsenal Sinus Tarsi Calc Plate, L	00812926029849
300-86-006	Arsenal Sinus Tarsi Calc Plate, R	00812926029856
300-86-007	Arsenal Sinus Support Plate	00812926029863
300-87-001	Arsenal Short Medial Column Plate, Left	00812926029870
300-87-002	Arsenal Short Medial Column Plate, Right	00812926029887
300-87-003	Arsenal Long Medial Column Plate, Left	00812926029894
300-87-004	Arsenal Long Medial Column Plate, Right	00812926029900
300-88-001	Arsenal L Plate, Left	00810028390270
300-88-002	Arsenal L Plate, Right	00810028390287
301-24-010	2.4mm x 10mm GL Screw	00812926022635
301-24-012	2.4mm x 12mm GL Screw	00812926022642
301-24-014	2.4mm x 14mm GL Screw	00812926022659
301-24-016	2.4mm x 16mm GL Screw	00812926022666
301-24-018	2.4mm x 18mm GL Screw	00812926022673
301-24-020	2.4mm x 20mm GL Screw	00812926022680
301-24-022	2.4mm x 22mm GL Screw	00812926022697
301-24-024	2.4mm x 24mm GL Screw	00812926022703
301-24-026	2.4mm x 26mm GL Screw	00812926022710
301-24-028	2.4mm x 28mm GL Screw	00812926022727
301-24-030	2.4mm x 30mm GL Screw	00812926022734
301-24-032	2.4mm x 32mm GL Screw	00812926022741
301-24-034	2.4mm x 34mm GL Screw	00812926022758
301-24-036	2.4mm x 36mm GL Screw	00812926022765
301-24-038	2.4mm x 38mm GL Screw	00812926022772
301-24-040	2.4mm x 40mm GL Screw	00812926022789
301-30-010	3.0mm x 10mm GL Screw	00812926022796
301-30-012	3.0mm x 12mm GL Screw	00812926022802
301-30-014	3.0mm x 14mm GL Screw	00812926022819
301-30-016	3.0mm x 16mm GL Screw	00812926022826
301-30-018	3.0mm x 18mm GL Screw	00812926022833
301-30-020	3.0mm x 20mm GL Screw	00812926022840
301-30-022	3.0mm x 22mm GL Screw	00812926022857
301-30-024	3.0mm x 24mm GL Screw	00812926022864
301-30-026	3.0mm x 26mm GL Screw	00812926022871
301-30-028	3.0mm x 28mm GL Screw	00812926022888

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
301-30-030	3.0mm x 30mm GL Screw	00812926022895
301-30-032	3.0mm x 32mm GL Screw	00812926022901
301-30-034	3.0mm x 34mm GL Screw	00812926022918
301-30-036	3.0mm x 36mm GL Screw	00812926022925
301-30-038	3.0mm x 38mm GL Screw	00812926022932
301-30-040	3.0mm x 40mm GL Screw	00812926022949
301-40-010	4.0mm x 10mm GL Screw	00812926022956
301-40-012	4.0mm x 12mm GL Screw	00812926022963
301-40-014	4.0mm x 14mm GL Screw	00812926022970
301-40-016	4.0mm x 16mm GL Screw	00812926022987
301-40-018	4.0mm x 18mm GL Screw	00812926022994
301-40-020	4.0mm x 20mm GL Screw	00812926023007
301-40-022	4.0mm x 22mm GL Screw	00812926023014
301-40-024	4.0mm x 24mm GL Screw	00812926023021
301-40-026	4.0mm x 26mm GL Screw	00812926023038
301-40-028	4.0mm x 28mm GL Screw	00812926023045
301-40-030	4.0mm x 30mm GL Screw	00812926023052
301-40-032	4.0mm x 32mm GL Screw	00812926023069
301-40-034	4.0mm x 34mm GL Screw	00812926023076
301-40-036	4.0mm x 36mm GL Screw	00812926023083
301-40-038	4.0mm x 38mm GL Screw	00812926023090
301-40-040	4.0mm x 40mm GL Screw	00812926023106
302-24-010	2.4mm x 10mm GL Locking Screw	00812926023113
302-24-012	2.4mm x 12mm GL Locking Screw	00812926023120
302-24-014	2.4mm x 14mm GL Locking Screw	00812926023137
302-24-016	2.4mm x 16mm GL Locking Screw	00812926023144
302-24-018	2.4mm x 18mm GL Locking Screw	00812926023151
302-24-020	2.4mm x 20mm GL Locking Screw	00812926023168
302-24-022	2.4mm x 22mm GL Locking Screw	00812926023175
302-24-024	2.4mm x 24mm GL Locking Screw	00812926023182
302-24-026	2.4mm x 26mm GL Locking Screw	00812926023199
302-24-028	2.4mm x 28mm GL Locking Screw	00812926023205
302-24-030	2.4mm x 30mm GL Locking Screw	00812926023212
302-24-032	2.4mm x 32mm GL Locking Screw	00812926023229
302-24-034	2.4mm x 34mm GL Locking Screw	00812926023236
302-24-036	2.4mm x 36mm GL Locking Screw	00812926023243
302-24-038	2.4mm x 38mm GL Locking Screw	00812926023250
302-24-040	2.4mm x 40mm GL Locking Screw	00812926023267
302-30-010	3.0mm x 10mm GL Locking Screw	00812926023274
302-30-012	3.0mm x 12mm GL Locking Screw	00812926023281
302-30-014	3.0mm x 14mm GL Locking Screw	00812926023298
302-30-016	3.0mm x 16mm GL Locking Screw	00812926023304
302-30-018	3.0mm x 18mm GL Locking Screw	00812926023311
302-30-020	3.0mm x 20mm GL Locking Screw	00812926023328
302-30-022	3.0mm x 22mm GL Locking Screw	00812926023335
302-30-024	3.0mm x 24mm GL Locking Screw	00812926023342

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
302-30-026	3.0mm x 26mm GL Locking Screw	00812926023359
302-30-028	3.0mm x 28mm GL Locking Screw	00812926023366
302-30-030	3.0mm x 30mm GL Locking Screw	00812926023373
302-30-032	3.0mm x 32mm GL Locking Screw	00812926023380
302-30-034	3.0mm x 34mm GL Locking Screw	00812926023397
302-30-036	3.0mm x 36mm GL Locking Screw	00812926023403
302-30-038	3.0mm x 38mm GL Locking Screw	00812926023410
302-30-040	3.0mm x 40mm GL Locking Screw	00812926023427
302-40-010	4.0mm x 10mm GL Locking Screw	00812926023434
302-40-012	4.0mm x 12mm GL Locking Screw	00812926023441
302-40-014	4.0mm x 14mm GL Locking Screw	00812926023458
302-40-016	4.0mm x 16mm GL Locking Screw	00812926023465
302-40-018	4.0mm x 18mm GL Locking Screw	00812926023472
302-40-020	4.0mm x 20mm GL Locking Screw	00812926023489
302-40-022	4.0mm x 22mm GL Locking Screw	00812926023496
302-40-024	4.0mm x 24mm GL Locking Screw	00812926023502
302-40-026	4.0mm x 26mm GL Locking Screw	00812926023519
302-40-028	4.0mm x 28mm GL Locking Screw	00812926023526
302-40-030	4.0mm x 30mm GL Locking Screw	00812926023533
302-40-032	4.0mm x 32mm GL Locking Screw	00812926023540
302-40-034	4.0mm x 34mm GL Locking Screw	00812926023557
302-40-036	4.0mm x 36mm GL Locking Screw	00812926023564
302-40-038	4.0mm x 38mm GL Locking Screw	00812926023571
302-40-040	4.0mm x 40mm GL Locking Screw	00812926023588
303-50-001	Smooth GL Screw Washer	00812926023595
303-50-002	Spiked GL Screw Washer	00812926023601
303-50-003	Large Smooth GL Screw Washer	00812926023618
303-50-004	12mm GL Spiked Washer	00812926023625
303-60-001	Smooth Gridlock Ankle Screw Washer	00812926026633
303-60-002	Spiked Gridlock Ankle Screw Washer	00812926026640
304-27-010	2.7mm x 10mm Gridlock Ankle Screw	00812926025513
304-27-012	2.7mm x 12mm Gridlock Ankle Screw	00812926025520
304-27-014	2.7mm x 14mm Gridlock Ankle Screw	00812926025537
304-27-016	2.7mm x 16mm Gridlock Ankle Screw	00812926025544
304-27-018	2.7mm x 18mm Gridlock Ankle Screw	00812926025551
304-27-020	2.7mm x 20mm Gridlock Ankle Screw	00812926025568
304-27-022	2.7mm x 22mm Gridlock Ankle Screw	00812926025575
304-27-024	2.7mm x 24mm Gridlock Ankle Screw	00812926025582
304-35-010	3.5mm x 10mm Gridlock Ankle Screw	00812926025599
304-35-012	3.5mm x 12mm Gridlock Ankle Screw	00812926025605
304-35-014	3.5mm x 14mm Gridlock Ankle Screw	00812926025612
304-35-016	3.5mm x 16mm Gridlock Ankle Screw	00812926025629
304-35-018	3.5mm x 18mm Gridlock Ankle Screw	00812926025636
304-35-020	3.5mm x 20mm Gridlock Ankle Screw	00812926025643
304-35-022	3.5mm x 22mm Gridlock Ankle Screw	00812926025650
304-35-024	3.5mm x 24mm Gridlock Ankle Screw	00812926025667



## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
304-35-026	3.5mm x 26mm Gridlock Ankle Screw	00812926025674
304-35-028	3.5mm x 28mm Gridlock Ankle Screw	00812926025681
304-35-030	3.5mm x 30mm Gridlock Ankle Screw	00812926025698
304-35-032	3.5mm x 32mm Gridlock Ankle Screw	00812926025704
304-35-034	3.5mm x 34mm Gridlock Ankle Screw	00812926025711
304-35-036	3.5mm x 36mm Gridlock Ankle Screw	00812926025728
304-35-038	3.5mm x 38mm Gridlock Ankle Screw	00812926025735
304-35-040	3.5mm x 40mm Gridlock Ankle Screw	00812926025742
304-35-042	3.5mm x 42mm Gridlock Ankle Screw	00812926025759
304-35-044	3.5mm x 44mm Gridlock Ankle Screw	00812926025766
304-35-046	3.5mm x 46mm Gridlock Ankle Screw	00812926025773
304-35-048	3.5mm x 48mm Gridlock Ankle Screw	00812926025780
304-35-050	3.5mm x 50mm Gridlock Ankle Screw	00812926025797
304-40-010	4.0mm x 10mm Gridlock Ankle Screw	00812926027692
304-40-012	4.0mm x 12mm Gridlock Ankle Screw	00812926025803
304-40-014	4.0mm x 14mm Gridlock Ankle Screw	00812926025810
304-40-016	4.0mm x 16mm Gridlock Ankle Screw	00812926025827
304-40-018	4.0mm x 18mm Gridlock Ankle Screw	00812926025834
304-40-020	4.0mm x 20mm Gridlock Ankle Screw	00812926025841
304-40-022	4.0mm x 22mm Gridlock Ankle Screw	00812926025858
304-40-024	4.0mm x 24mm Gridlock Ankle Screw	00812926025865
304-40-026	4.0mm x 26mm Gridlock Ankle Screw	00812926025872
304-40-028	4.0mm x 28mm Gridlock Ankle Screw	00812926025889
304-40-030	4.0mm x 30mm Gridlock Ankle Screw	00812926025896
304-40-032	4.0mm x 32mm Gridlock Ankle Screw	00812926025902
304-40-034	4.0mm x 34mm Gridlock Ankle Screw	00812926025919
304-40-036	4.0mm x 36mm Gridlock Ankle Screw	00812926025926
304-40-038	4.0mm x 38mm Gridlock Ankle Screw	00812926025933
304-40-040	4.0mm x 40mm Gridlock Ankle Screw	00812926025940
304-40-042	4.0mm x 42mm Gridlock Ankle Screw	00812926025957
304-40-044	4.0mm x 44mm Gridlock Ankle Screw	00812926025964
304-40-046	4.0mm x 46mm Gridlock Ankle Screw	00812926025971
304-40-048	4.0mm x 48mm Gridlock Ankle Screw	00812926025988
304-40-050	4.0mm x 50mm Gridlock Ankle Screw	00812926025995
305-27-010	2.7mm x 10mm Gridlock Ankle Locking Screw	00812926026008
305-27-012	2.7mm x 12mm Gridlock Ankle Locking Screw	00812926026015
305-27-014	2.7mm x 14mm Gridlock Ankle Locking Screw	00812926026022
305-27-016	2.7mm x 16mm Gridlock Ankle Locking Screw	00812926026039
305-27-018	2.7mm x 18mm Gridlock Ankle Locking Screw	00812926026046
305-27-020	2.7mm x 20mm Gridlock Ankle Locking Screw	00812926026053
305-27-022	2.7mm x 22mm Gridlock Ankle Locking Screw	00812926026060
305-27-024	2.7mm x 24mm Gridlock Ankle Locking Screw	00812926026077
305-35-010	3.5mm x 10mm Gridlock Ankle Locking Screw	00812926026084
305-35-012	3.5mm x 12mm Gridlock Ankle Locking Screw	00812926026091
305-35-014	3.5mm x 14mm Gridlock Ankle Locking Screw	00812926026107
305-35-016	3.5mm x 16mm Gridlock Ankle Locking Screw	00812926026114

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
305-35-018	3.5mm x 18mm Gridlock Ankle Locking Screw	00812926026121
305-35-020	3.5mm x 20mm Gridlock Ankle Locking Screw	00812926026138
305-35-022	3.5mm x 22mm Gridlock Ankle Locking Screw	00812926026145
305-35-024	3.5mm x 24mm Gridlock Ankle Locking Screw	00812926026152
305-35-026	3.5mm x 26mm Gridlock Ankle Locking Screw	00812926026169
305-35-028	3.5mm x 28mm Gridlock Ankle Locking Screw	00812926026176
305-35-030	3.5mm x 30mm Gridlock Ankle Locking Screw	00812926026183
305-35-032	3.5mm x 32mm Gridlock Ankle Locking Screw	00812926026190
305-35-034	3.5mm x 34mm Gridlock Ankle Locking Screw	00812926026206
305-35-036	3.5mm x 36mm Gridlock Ankle Locking Screw	00812926026213
305-35-038	3.5mm x 38mm Gridlock Ankle Locking Screw	00812926026220
305-35-040	3.5mm x 40mm Gridlock Ankle Locking Screw	00812926026237
305-35-042	3.5mm x 42mm Gridlock Ankle Locking Screw	00812926026244
305-35-044	3.5mm x 44mm Gridlock Ankle Locking Screw	00812926026251
305-35-046	3.5mm x 46mm Gridlock Ankle Locking Screw	00812926026268
305-35-048	3.5mm x 48mm Gridlock Ankle Locking Screw	00812926026275
305-35-050	3.5mm x 50mm Gridlock Ankle Locking Screw	00812926026282
305-40-010	4.0mm x 10mm Gridlock Ankle Locking Screw	00812926026299
305-40-012	4.0mm x 12mm Gridlock Ankle Locking Screw	00812926026305
305-40-014	4.0mm x 14mm Gridlock Ankle Locking Screw	00812926026312
305-40-016	4.0mm x 16mm Gridlock Ankle Locking Screw	00812926026329
305-40-018	4.0mm x 18mm Gridlock Ankle Locking Screw	00812926026336
305-40-020	4.0mm x 20mm Gridlock Ankle Locking Screw	00812926026343
305-40-022	4.0mm x 22mm Gridlock Ankle Locking Screw	00812926026350
305-40-024	4.0mm x 24mm Gridlock Ankle Locking Screw	00812926026367
305-40-026	4.0mm x 26mm Gridlock Ankle Locking Screw	00812926026374
305-40-028	4.0mm x 28mm Gridlock Ankle Locking Screw	00812926026381
305-40-030	4.0mm x 30mm Gridlock Ankle Locking Screw	00812926026398
305-40-032	4.0mm x 32mm Gridlock Ankle Locking Screw	00812926026404
305-40-034	4.0mm x 34mm Gridlock Ankle Locking Screw	00812926026411
305-40-036	4.0mm x 36mm Gridlock Ankle Locking Screw	00812926026428
305-40-038	4.0mm x 38mm Gridlock Ankle Locking Screw	00812926026435
305-40-040	4.0mm x 40mm Gridlock Ankle Locking Screw	00812926026442
305-40-042	4.0mm x 42mm Gridlock Ankle Locking Screw	00812926026459
305-40-044	4.0mm x 44mm Gridlock Ankle Locking Screw	00812926026466
305-40-046	4.0mm x 46mm Gridlock Ankle Locking Screw	00812926026473
305-40-048	4.0mm x 48mm Gridlock Ankle Locking Screw	00812926026480
305-40-050	4.0mm x 50mm Gridlock Ankle Locking Screw	00812926026497
306-40-040	4.0mm x 40mm Gridlock Ankle Syndesmotoc Screw	00812926026503
306-40-042	4.0mm x 42mm Gridlock Ankle Syndesmotoc Screw	00812926026510
306-40-044	4.0mm x 44mm Gridlock Ankle Syndesmotoc Screw	00812926026527
306-40-046	4.0mm x 46mm Gridlock Ankle Syndesmotoc Screw	00812926026534
306-40-048	4.0mm x 48mm Gridlock Ankle Syndesmotoc Screw	00812926026541
306-40-050	4.0mm x 50mm Gridlock Ankle Syndesmotoc Screw	00812926026558
306-40-052	4.0mm x 52mm Gridlock Ankle Syndesmotoc Screw	00812926026565
306-40-054	4.0mm x 54mm Gridlock Ankle Syndesmotoc Screw	00812926026572

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
306-40-056	4.0mm x 56mm Gridlock Ankle Syndesmotic Screw	00812926026589
306-40-058	4.0mm x 58mm Gridlock Ankle Syndesmotic Screw	00812926026596
306-40-060	4.0mm x 60mm Gridlock Ankle Syndesmotic Screw	00812926026602
306-40-062	4.0mm x 62mm Gridlock Ankle Syndesmotic Screw	00812926026619
306-40-064	4.0mm x 64mm Gridlock Ankle Syndesmotic Screw	00812926026626
306-41-040	4.0mm x 40mm Gridlock Ankle Syndesmotic Screw LP	00812926027234
306-41-042	4.0mm x 42mm Gridlock Ankle Syndesmotic Screw LP	00812926027241
306-41-044	4.0mm x 44mm Gridlock Ankle Syndesmotic Screw LP	00812926027258
306-41-046	4.0mm x 46mm Gridlock Ankle Syndesmotic Screw LP	00812926027265
306-41-048	4.0mm x 48mm Gridlock Ankle Syndesmotic Screw LP	00812926027272
306-41-050	4.0mm x 50mm Gridlock Ankle Syndesmotic Screw LP	00812926027289
306-41-052	4.0mm x 52mm Gridlock Ankle Syndesmotic Screw LP	00812926027296
306-41-054	4.0mm x 54mm Gridlock Ankle Syndesmotic Screw LP	00812926027302
306-41-056	4.0mm x 56mm Gridlock Ankle Syndesmotic Screw LP	00812926027319
306-41-058	4.0mm x 58mm Gridlock Ankle Syndesmotic Screw LP	00812926027326
306-41-060	4.0mm x 60mm Gridlock Ankle Syndesmotic Screw LP	00812926027333
306-41-062	4.0mm x 62mm Gridlock Ankle Syndesmotic Screw LP	00812926027340
306-41-064	4.0mm x 64mm Gridlock Ankle Syndesmotic Screw LP	00812926027357
307-22-006	2.2mm x 6mm Arsenal Screw	00810028390294
307-22-008	2.2mm x 8mm Arsenal Screw	00810028390300
307-22-010	2.2mm x 10mm Arsenal Screw	00810028390317
307-22-012	2.2mm x 12mm Arsenal Screw	00810028390324
307-22-014	2.2mm x 14mm Arsenal Screw	00810028390331
307-22-016	2.2mm x 16mm Arsenal Screw	00810028390348
307-22-018	2.2mm x 18mm Arsenal Screw	00810028390355
307-22-020	2.2mm x 20mm Arsenal Screw	00810028390362
307-22-022	2.2mm x 22mm Arsenal Screw	00810028390379
307-22-024	2.2mm x 24mm Arsenal Screw	00810028390386
307-22-026	2.2mm x 26mm Arsenal Screw	00810028390393
307-22-028	2.2mm x 28mm Arsenal Screw	00810028390409
307-22-030	2.2mm x 30mm Arsenal Screw	00810028390416
307-27-008	2.7mm x 8mm Arsenal Screw	00810028390423
307-27-010	2.7mm x 10mm Arsenal Screw	00810028390430
307-27-012	2.7mm x 12mm Arsenal Screw	00810028390447
307-27-014	2.7mm x 14mm Arsenal Screw	00810028390454
307-27-016	2.7mm x 16mm Arsenal Screw	00810028390461
307-27-018	2.7mm x 18mm Arsenal Screw	00810028390478
307-27-020	2.7mm x 20mm Arsenal Screw	00810028390485
307-27-022	2.7mm x 22mm Arsenal Screw	00810028390492
307-27-024	2.7mm x 24mm Arsenal Screw	00810028390508
307-27-026	2.7mm x 26mm Arsenal Screw	00810028390515
307-27-028	2.7mm x 28mm Arsenal Screw	00810028390522
307-27-030	2.7mm x 30mm Arsenal Screw	00810028390539
307-27-032	2.7mm x 32mm Arsenal Screw	00810028390546
307-27-034	2.7mm x 34mm Arsenal Screw	00810028390553
307-27-036	2.7mm x 36mm Arsenal Screw	00810028390560

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
307-27-038	2.7mm x 38mm Arsenal Screw	00810028390577
307-27-040	2.7mm x 40mm Arsenal Screw	00810028390584
307-35-008	3.5mm x 8mm Arsenal Screw	00810028390591
307-35-010	3.5mm x 10mm Arsenal Screw	00810028390607
307-35-012	3.5mm x 12mm Arsenal Screw	00810028390614
307-35-014	3.5mm x 14mm Arsenal Screw	00810028390621
307-35-016	3.5mm x 16mm Arsenal Screw	00810028390638
307-35-018	3.5mm x 18mm Arsenal Screw	00810028390645
307-35-020	3.5mm x 20mm Arsenal Screw	00810028390652
307-35-022	3.5mm x 22mm Arsenal Screw	00810028390669
307-35-024	3.5mm x 24mm Arsenal Screw	00810028390676
307-35-026	3.5mm x 26mm Arsenal Screw	00810028390683
307-35-028	3.5mm x 28mm Arsenal Screw	00810028390690
307-35-030	3.5mm x 30mm Arsenal Screw	00810028390706
307-35-032	3.5mm x 32mm Arsenal Screw	00810028390713
307-35-034	3.5mm x 34mm Arsenal Screw	00810028390720
307-35-036	3.5mm x 36mm Arsenal Screw	00810028390737
307-35-038	3.5mm x 38mm Arsenal Screw	00810028390744
307-35-040	3.5mm x 40mm Arsenal Screw	00810028390751
307-35-042	3.5mm x 42mm Arsenal Screw	00810028390768
307-35-044	3.5mm x 44mm Arsenal Screw	00810028390775
307-35-046	3.5mm x 46mm Arsenal Screw	00810028390782
307-35-048	3.5mm x 48mm Arsenal Screw	00810028390799
307-35-050	3.5mm x 50mm Arsenal Screw	00810028390805
308-22-006	2.2mm x 6mm Arsenal Locking Screw	00810028390812
308-22-008	2.2mm x 8mm Arsenal Locking Screw	00810028390829
308-22-010	2.2mm x 10mm Arsenal Locking Screw	00810028390836
308-22-012	2.2mm x 12mm Arsenal Locking Screw	00810028390843
308-22-014	2.2mm x 14mm Arsenal Locking Screw	00810028390850
308-22-016	2.2mm x 16mm Arsenal Locking Screw	00810028390867
308-22-018	2.2mm x 18mm Arsenal Locking Screw	00810028390874
308-22-020	2.2mm x 20mm Arsenal Locking Screw	00810028390881
308-22-022	2.2mm x 22mm Arsenal Locking Screw	00810028390898
308-22-024	2.2mm x 24mm Arsenal Locking Screw	00810028390904
308-22-026	2.2mm x 26mm Arsenal Locking Screw	00810028390911
308-22-028	2.2mm x 28mm Arsenal Locking Screw	00810028390928
308-22-030	2.2mm x 30mm Arsenal Locking Screw	00810028390935
308-27-008	2.7mm x 8mm Arsenal Locking Screw	00810028390942
308-27-010	2.7mm x 10mm Arsenal Locking Screw	00810028390959
308-27-012	2.7mm x 12mm Arsenal Locking Screw	00810028390966
308-27-014	2.7mm x 14mm Arsenal Locking Screw	00810028390973
308-27-016	2.7mm x 16mm Arsenal Locking Screw	00810028390980
308-27-018	2.7mm x 18mm Arsenal Locking Screw	00810028390997
308-27-020	2.7mm x 20mm Arsenal Locking Screw	00810028391000
308-27-022	2.7mm x 22mm Arsenal Locking Screw	00810028391017
308-27-024	2.7mm x 24mm Arsenal Locking Screw	00810028391024

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
308-27-026	2.7mm x 26mm Arsenal Locking Screw	00810028391031
308-27-028	2.7mm x 28mm Arsenal Locking Screw	00810028391048
308-27-030	2.7mm x 30mm Arsenal Locking Screw	00810028391055
308-27-032	2.7mm x 32mm Arsenal Locking Screw	00810028391062
308-27-034	2.7mm x 34mm Arsenal Locking Screw	00810028391079
308-27-036	2.7mm x 36mm Arsenal Locking Screw	00810028391086
308-27-038	2.7mm x 38mm Arsenal Locking Screw	00810028391093
308-27-040	2.7mm x 40mm Arsenal Locking Screw	00810028391109
308-35-008	3.5mm x 8mm Arsenal Locking Screw	00810028391116
308-35-010	3.5mm x 10mm Arsenal Locking Screw	00810028391123
308-35-012	3.5mm x 12mm Arsenal Locking Screw	00810028391130
308-35-014	3.5mm x 14mm Arsenal Locking Screw	00810028391147
308-35-016	3.5mm x 16mm Arsenal Locking Screw	00810028391154
308-35-018	3.5mm x 18mm Arsenal Locking Screw	00810028391161
308-35-020	3.5mm x 20mm Arsenal Locking Screw	00810028391178
308-35-022	3.5mm x 22mm Arsenal Locking Screw	00810028391185
308-35-024	3.5mm x 24mm Arsenal Locking Screw	00810028391192
308-35-026	3.5mm x 26mm Arsenal Locking Screw	00810028391208
308-35-028	3.5mm x 28mm Arsenal Locking Screw	00810028391215
308-35-030	3.5mm x 30mm Arsenal Locking Screw	00810028391222
308-35-032	3.5mm x 32mm Arsenal Locking Screw	00810028391239
308-35-034	3.5mm x 34mm Arsenal Locking Screw	00810028391246
308-35-036	3.5mm x 36mm Arsenal Locking Screw	00810028391253
308-35-038	3.5mm x 38mm Arsenal Locking Screw	00810028391260
308-35-040	3.5mm x 40mm Arsenal Locking Screw	00810028391277
308-35-042	3.5mm x 42mm Arsenal Locking Screw	00810028391284
308-35-044	3.5mm x 44mm Arsenal Locking Screw	00810028391291
308-35-046	3.5mm x 46mm Arsenal Locking Screw	00810028391307
308-35-048	3.5mm x 48mm Arsenal Locking Screw	00810028391314
308-35-050	3.5mm x 50mm Arsenal Locking Screw	00810028391321
310-00-001	Large Bone Clamp	00812926023632
310-00-002	Small Bone Clamp	00812926023649
310-00-003	5" Joint Distraction Plier	00812926023656
310-00-004	4" Osteotomy Distraction Plier	00812926023663
310-00-006	Plate Bending Irons	00812926023670
310-00-008	Large Bending Pliers	00812926023687
310-00-009	GL Calibrated Screw Sleeve	00812926028040
310-00-010	Plate Screw Depth Gauge	00812926023694
310-00-019	GL Plate Forceps	00812926024066
310-00-020	6" Distraction Pliers	00812926028019
310-00-021	8" Distraction Pliers	00812926028026
310-01-001	Universal IM Guide Pin	00812926023755
310-01-002	Cannulated Reamer Guide	00812926023762
310-01-016	16mm Concave Reamer Head	00812926023779
310-01-018	18mm Concave Reamer Head	00812926023793
310-01-020	20mm Concave Reamer Head	00812926023816

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
310-01-022	22mm Concave Reamer Head	00812926023830
310-01-100	Gridlock Sterilization Tray Assembly	00812926023700
310-01-101	Gridlock Plating System Screw Module Base Assembly	00812926023717
310-01-102	Gridlock VLC MPX L Plating Module Base Assembly	00812926023724
310-01-103	Gridlock Instrument Module Base Assembly	00812926023731
310-01-104	Gridlock Specialty Plate Module Base Assembly	00812926023748
310-02-016	16mm Convex Reamer Head	00812926023786
310-02-018	18mm Convex Reamer Head	00812926023809
310-02-020	20mm Convex Reamer Head	00812926023823
310-02-022	22mm Convex Reamer Head	00812926023847
310-24-001	2.4mm GL Screw Drill Bit	00812926023854
310-24-005	2.4mm Variable/Static Angle Drill Guide	00812926023984
310-24-006	2.4mm Drill Sleeves	00812926024011
310-24-007	2.4mm Compression/Neutral Drill Guide	00812926024042
310-24-008	2.4mm GL Static Threaded Drill Guide	00812926025414
310-24-009	2.4mm Variable Angle Threaded Drill Guide	00812926028293
310-25-001	2.4mm GL Screw Overdrill Bit	00812926023885
310-30-001	3.0mm GL Screw Drill Bit	00812926023861
310-30-002	2.4/3.0mm GL Screw Countersink	00812926023915
310-30-003	2.4/3.0/4.0mm GL Screw Driver Bit	00812926023977
310-30-005	3.0mm Variable/Static Angle Drill Guide	00812926023991
310-30-006	3.0mm Drill Sleeves	00812926024028
310-30-007	3.0mm Compression/Neutral Drill Guide	00812926024059
310-30-008	3.0mm GL Static Threaded Drill Guide	00812926025421
310-30-009	3.0mm Variable Angle Threaded Drill Guide	00812926028309
310-31-001	3.0mm GL Screw Overdrill Bit	00812926023892
310-31-003	Gridlock Driver Bit, Short	00812926027548
310-40-001	4.0mm GL Pilot Drill Bit	00812926023878
310-40-002	2.4/3.0/4.0mm GL Screw Countersink	00812926023922
310-40-005	4.0mm Variable/Static Angle Drill Guide	00812926024004
310-40-006	4.0mm Drill Sleeves	00812926024035
310-40-009	4.0mm Variable Angle Threaded Drill Guide	00812926028316
310-41-001	4.0mm GL Overdrill Bit	00812926023908
310-60-005	0.062" x 4.5" Olive Wire Partially Threaded	00812926023939
310-99-033	Gridlock Plating System Screw Module Lid	00812926027920
310-99-034	Gridlock Plating System Instrument Module Lid	00812926027937
310-99-105	Gridlock VLC MPX L Plating Module Lid	00812926027944
310-99-107	Gridlock Specialty Plate Module Lid	00812926027951
311-00-001	Small Joint Reamer Sterilization Tray Base Assembly	00812926027463
311-01-016	16mm Concave Reamer	00812926027364
311-01-018	18mm Concave Reamer	00812926027388
311-01-020	20mm Concave Reamer	00812926027401
311-01-022	22mm Concave Reamer	00812926027425
311-01-024	24mm Concave Reamer	00812926027449
311-02-016	16mm Convex Reamer	00812926027371
311-02-018	18mm Convex Reamer	00812926027395

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
311-02-020	20mm Convex Reamer	00812926027418
311-02-022	22mm Convex Reamer	00812926027432
311-02-024	24mm Convex Reamer	00812926027456
311-99-011	Small Joint Reaming System Lid	00812926027968
317-00-001	Minimally Invasive Bunion Plate System Sterilization Tray	00812926027722
317-00-002	Plate Placement Guide, Right	00812926027746
317-00-003	Plate Placement Guide, Left	00812926027753
317-00-004	MIB Threaded Drill Guide	00812926027760
317-00-005	MIB Interfrag Drill Guide	00812926027777
317-00-006	Broach	00812926027784
317-00-007	Tissue Retractor	00812926027791
317-24-001	2.4mm Calibrated Drill Bit	00812926027807
317-99-011	Minimally Invasive Bunion Plate System Sterilization Tray Lid	00812926027739
320-00-001	Large Serrated Bone Clamp	00812926026695
320-00-003	Syndesmotic Bone Clamp	00812926026701
320-00-005	Fixed Driver Handle Gridlock Ankle Plating System	00812926027005
320-00-006	Gridlock Ankle Plate Bending Irons	00812926026718
320-00-009	GL Calibrated Screw Sleeve, Ankle	00812926028057
320-00-010	Gridlock Ankle Plating Depth Gauge	00812926026725
320-00-011	8mm Hohmann Retractor	00812926026732
320-00-012	15mm Hohmann Retractor	00812926026749
320-00-013	Bone Pick	00812926026756
320-00-014	Periosteal Elevator	00812926026763
320-00-016	GL Ankle Screw Remover	00812926026770
320-00-017	Hook Drill Guide	00812926027654
320-00-019	Guide Handle	00812926028798
320-00-020	Hook Tamp	00812926028804
320-01-100	Gridlock Ankle Plate System Tray Assembly	00812926026657
320-01-101	Ankle Screw Module Total Base Assembly	00812926026664
320-01-102	Gridlock Ankle Plate Module Base Assembly	00812926026671
320-01-103	Ankle Instrument Module Base Assembly	00812926026688
320-01-104	Ankle Plate Module, Ext.	00812926028286
320-01-999	Ankle Plate Module Insert	00812926027975
320-27-001	2.7mm Gridlock Ankle Screw Drill Bit	00812926026794
320-27-005	2.7mm Ankle Variable/Static Angle Drill Guide	00812926026909
320-27-006	2.7mm Ankle Drill Sleeves	00812926026947
320-27-008	2.7mm Ankle Threaded Static Drill Guide	00812926026978
320-27-009	2.7mm Variable Angle Threaded Drill Guide	00812926028323
320-27-011	2.7mm Gridlock Ankle Screw Overdrill Bit	00812926026824
320-35-001	3.5mm Gridlock Ankle Screw Drill Bit	00812926026800
320-35-002	Gridlock Ankle Screw Countersink	00812926026862
320-35-003	Gridlock Ankle Screw Driver Bit	00812926026893
320-35-005	3.5mm Ankle Variable/Static Angle Drill Guide	00812926026916
320-35-006	3.5mm Ankle Drill Sleeves	00812926026954
320-35-008	3.5mm Ankle Threaded Static Drill Guide	00812926026985
320-35-009	3.5mm Variable Angle Threaded Drill Guide	00812926028330

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
320-35-011	3.5mm Gridlock Ankle Screw Overdrill Bit	00812926026831
320-36-003	Gridlock Ankle Driver Bit, Short	00812926027555
320-40-001	4.0mm Gridlock Ankle Screw Drill Bit	00812926026817
320-40-002	Gridlock Ankle Washer Countersink	00812926026879
320-40-005	4.0mm Ankle Variable/Static Angle Drill Guide	00812926026923
320-40-006	4.0mm Ankle Drill Sleeves	00812926026961
320-40-008	4.0mm Ankle Threaded Static Drill Guide	00812926026992
320-40-009	4.0mm Variable Angle Threaded Drill Guide	00812926028347
320-40-011	4.0mm Gridlock Ankle Screw Overdrill Bit	00812926026848
320-40-012	Syndesmotic Screw Drill Bit	00812926026855
320-40-015	Syndesmotic Screw Drill Guide	00812926026930
320-40-016	Syndesmotic Screw Removal Bit	00812926026787
320-60-005	0.062" x 4.65" GL Olive Wire Partially Threaded	00812926026886
320-99-111	Gridlock Ankle Screw Module Lid	00812926027982
320-99-112	Gridlock Ankle Plate Module Lid	00812926027999
320-99-113	Gridlock Ankle Instrument Module Lid	00812926028002
330-00-001	Calibrated Spreader	00812926029528
330-00-002	Small Bone Distractor	00812926029542
330-00-003	Small Bone Compressor	00812926029559
330-00-004	McGlamry Elevator	00812926029566
330-00-005	Arsenal Variable Angle Drill Guide	00812926029740
330-00-006	Arsenal Compression Drill Guide	00812926029757
330-00-007	Arsenal Static Angle Drill Guide	00812926029764
330-00-008	2.7mm/3.5mm Arsenal Threaded Static Drill Guide	00812926029795
330-00-009	Arsenal Screw Remover	00810028391338
330-00-010	Large Bending Pliers with Springs	00810028391352
330-01-100	Arsenal Plating System Tray	00812926029801
330-07-003	2.2mm Arsenal Screw Driver Bit	00812926029658
330-10-003	2.7mm/3.5mm Arsenal Screw Driver Bit	00812926029665
330-22-001	2.2mm Arsenal Screw Drill Bit	00812926029672
330-22-005	2.2mm Arsenal Screw Depth Gauge	00812926029597
330-22-006	2.2mm Arsenal Drill Sleeves	00812926029733
330-22-011	2.2mm Arsenal Screw Overdrill Bit	00812926029702
330-24-004	0.094" x 6" K-Wire Standard	00812926029610
330-27-001	2.7mm Arsenal Screw Drill Bit	00812926029689
330-27-002	2.7mm/3.5mm Arsenal Screw Countersink	00812926029641
330-27-005	2.7mm/3.5mm Arsenal Screw Depth Gauge	00812926029603
330-27-006	2.7mm Arsenal Drill Sleeves	00812926029771
330-27-011	2.7mm Arsenal Screw Overdrill Bit	00812926029719
330-35-001	3.5mm Arsenal Screw Drill Bit	00812926029696
330-35-006	3.5mm Arsenal Drill Sleeves	00812926029788
330-35-011	3.5mm Arsenal Screw Overdrill Bit	00812926029726
330-60-005	0.062" x 3.25" Olive Wire Partially Threaded	00812926029627
330-86-001	Calc Slide Plate Placement Guide	00812926029634
500-08-001	8mm Sterile Staple Kit Assembly	00812926027500
500-08-101	8mm x 8mm x 8mm Staple	00812926028217



## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)
500-10-001	10mm Sterile Staple Kit Assembly	00812926027517
500-10-101	10mm x 10mm x 10mm Staple	00812926028224
500-12-101	12mm x 12mm x 12mm Staple	00812926028606
500-15-101	15mm x 15mm x 15mm Staple	00812926028231
500-15-102	15mm x 12mm x 12mm Staple	00812926028248
500-18-101	18mm x 18mm x 18mm Staple	00812926028255
500-18-102	18mm x 15mm x 15mm Staple	00812926028262
500-20-101	20mm x 20mm x 20mm Staple	00812926028279
500-20-201	20mm x 20mm x 20mm High Strength Staple	00812926028613
500-25-101	25mm x 20mm x 20mm Staple	00812926028620
510-00-001	Sniper Staple System Sterilization Tray Base Assembly	00812926028194
510-00-003	8mm, 10mm Staple Spreader	00812926028149
510-00-004	15mm, 18mm, 20mm Staple Spreader	00812926028156
510-00-005	Staple Tamp	00812926028187
510-00-006	8-12mm Staple Spreader	00812926028637
510-00-007	15-25mm Staple Spreader	00812926028644
510-00-010	Lever Assist	00812926027845
510-08-101	8mm Drill Guide Assembly	00812926028064
510-10-002	8mm/10mm Staple Pilot Drill	00812926028118
510-10-101	10mm Drill Guide Assembly	00812926028071
510-12-001	12mm Drill Guide Assembly	00812926028668
510-15-001	15mm Drill Guide Assembly	00812926028088
510-16-001	1.6mm Post / Pilot Drill	00812926028521
510-16-002	1.6mm Pilot Drill	00812926028538
510-16-005	1.6mm Locating Pin	00812926028545
510-18-001	18mm Drill Guide Assembly	00812926028095
510-20-001	20mm Drill Guide Assembly	00812926028101
510-20-002	2.0mm Calibrated Pilot Drill	00812926028125
510-20-005	2.0mm Locating Pin	00812926028163
510-20-101	20mm, HS, Drill Guide Assembly	00812926028675
510-25-001	25mm Drill Guide Assembly	00812926028682
510-25-002	2.5mm Calibrated Pilot Drill	00812926028132
510-25-005	2.5mm Locating Pin	00812926028170
510-30-002	3.0mm Calibrated Pilot Drill	00812926028651
510-30-005	3.0mm Locating Pin	00812926028590
510-99-101	Sniper Staple System Sterilization Tray Lid	00812926028200